http://www.lenstec.com  
ISO13485 Registered Company

LENSTEC (BARBADOS) INC., Airport Commercial Centre, Pilgrim Road, Christ Church BB1/092, Barbados

+1 (246) 420-6795 +1 (246) 420-6797 lenstecbarbados@lenstec.com

EC REP CMI spol. s.r.o. Trenčianska 47,821 09 Bratislava, Slovakia

+421 (2) 52631441 +421 (2) 52922882 lenstec@cmi.sk

UK SALES: LENSTEC (BARBADOS) INC., Lenstec House, Unit 8, Mariner Court, Calder Park, Wakefield, WF4 3FL, UK



+44 (0) 1924 382 678 +44 (0) 1924 850 454 lenstecuk@lenstec.com

WORLD HEADQUARTERS: LENSTEC, INC., 1765 Commerce Ave N, Saint Petersburg, Florida 33716, U.S.A.


+1 (727) 571-2272 +1 (727) 571-1792 lenstec@lenstec.com

INSTRUCTIONS FOR USE • MODE D'EMPLOI • INSTRUCCIONES DE USO • ISTRUZIONI PER L'USO • GEBRAUCHSANWEISUNG • GEBRUIKSAANWIJZING • INSTRUÇÕES DE UTILIZAÇÃO  
• NÁVOD K POUŽITÍ • POKYNY K POUŽITIU • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • KULLANIM TALÍMATLARI  
• إرشادات الاستخدام • 使用説明

## MASTER

PCIOI

ACIOI

DocuSigned by:

LB1 Document Admin

94B46562700E486

11-Nov-2021

DCR5906/  
CRP211

EN	INSTRUCTIONS FOR USE	2
FR	MODE D'EMPLOI	3
ES	INSTRUCCIONES DE USO	5
IT	ISTRUZIONI PER L'USO	7
DE	GEBRAUCHSANWEISUNG	9
PT	INSTRUÇÕES DE UTILIZAÇÃO	11
CS	NÁVOD K POUŽITÍ	12
SK	POKYNY K POUŽITIU	14
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	15
TR	KULLANIM TALÍMATLARI	18
AR	إرشادات الاستخدام	20
ZH-Hans	使用説明	22

CE  
2797

STERILE EO

Rx only

PI12 Rev 13

**IMPORTANT NOTICE**

It is highly recommended that the surgeon adheres to the recommendations, contraindications and warnings outlined in these instructions.

**INTRODUCTION**

Lenstec PMMA intraocular lenses are optical implants used in the visual correction of cataracts. They are single-piece, equiconvex, ultraviolet absorbing intraocular lenses with step-vaulted, modified C loop or plate haptics. The PMMA IOLs are manufactured completely from Perspex CQ-UV medical grade polymethylmethacrylate (PMMA). This includes IOLs intended for placement in both the anterior and posterior chambers.

**HOW THE DEVICE IS USED****INDICATIONS**

Lenstec PMMA intraocular lenses are optical implants used in the replacement of cataractous natural lenses, and are intended to be placed in either the anterior or posterior chamber, based on model.

**RISK ANALYSIS**

The potential risks of implanting the Lenstec LS-Series IOLs are as follows:

Endophthalmitis, Toxic anterior segment syndrome (TASS), Anterior capsule fibrosis, Uveitis glaucoma hyphaema syndrome, Iritis, Iris capture, Cystoid macular edema, Corneal stromal edema, Posterior capsular contraction & lens deformation, Capsular damage, Decentration/tilt (small optic), Elevated IOP

The potential risks of implanting the Lenstec LA-Series IOLs are as follows:

Corneal decompensation, Pupil distortion, Peripheral anterior synechiae, Glaucoma, Uveitis-glaucoma-hyphaema syndrome, Endophthalmitis, Toxic anterior segment syndrome (TASS), Anterior capsule fibrosis, Iritis, Iris capture, Cystoid macular edema, Corneal stromal edema, Capsular damage, Decentration/tilt (small optic), Elevated IOP

Individuals in the following health categories may be at a higher risk of experiencing secondary IOL calcification (surface deposits on the lens):

- Diabetes
  - ◆ Associated:
    - \* Retinal detachment
    - \* Vitreous detachment
    - \* Vitrectomy
    - \* Diabetic retinopathy
    - \* Diabetic maculopathy
- Glaucoma
- Underwent DMEK or DSEK (multiple procedures further increase the risk)
- Hypertension
- High cholesterol

**MASTER****CONTRAINDICATIONS**

Outside of general contraindications for ocular surgery, the following specific contraindications apply:

Uncontrolled glaucoma, microphthalmia, chronic severe uveitis, retinal detachment, corneal decompensation, diabetic retinopathy, iris atrophy, perioperative complications, potentially foreseeable post-operative complications and other conditions which an ophthalmic surgeon might identify based on their experience.

**WARNINGS**

As with any surgical procedure, there is potential risk involved. Potential complications accompanying intraocular lens implantation may include, but are not limited to the following:

Lens Dislocation, Vitreous Loss, Corneal Endothelial Damage, Pupillary Block, Non-pigment Precipitates, Secondary Glaucoma, Cystoid Macular Edema, Iris Prolapse, Infection, Vitreous Wick Syndrome, Retinal Detachment, Papillary Membrane.

Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

- Recurrent severe anterior or posterior segment inflammation or uveitis.
- Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment diseases.
- Surgical difficulties at the time of cataract extraction which might increase the potential for complications (e.g. persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
- A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
- Circumstances that would result in damage to the endothelium during implantation.
- Suspected microbial infection.
- Children under the age of 2 years are not suitable candidates for intraocular lenses.
- Re-use of the IOL is strictly prohibited, as it raises serious safety and efficacy concerns.
  - LENSTEC does not provide cleaning/sterilization instructions. An improperly cleaned and/or sterilized IOL can cause significant damage to a patient's vision, due in part to cross contamination induced infection.
  - Once removed from its original packaging, the IOL can lose traceability. In the event an IOL is re-used, it is unlikely the user will know the correct expiry date, serial number or dioptric power.
  - LENSTEC can not guarantee stability or proper function of either haptic or optic portions in the event that an IOL is re-used. Failure of either of these components can render the IOL ineffective.
- The effectiveness of this ultraviolet absorbing lens in reducing the incidence of retinal disorders has not been established.
- Patients with ocular pathology may not achieve the visual acuity and/or have increased complications compared to patients without such pathology. Physicians should explore the use of alternative methods of aphakia correction in these patients and should consider lens implantation only if alternative treatments are deemed unsatisfactory to meet the needs of the patient.
- Secondary glaucoma has been reported occasionally in patients with preexisting glaucoma who received lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored postoperatively.
- Hyphema, secondary glaucoma, papillary block, cystic membrane formation and vitritis have been reported at increased rates in

patients who have surgical complications associated with cataract extraction procedure.

- Patients who have operative complications should be monitored postoperatively for the occurrence of these complications.

## PRECAUTIONS

- Do not use if the packaging has been opened or damaged.
- Do not reuse the lens.
- Do not sterilize. Lenses requiring re-sterilization should be returned to LENSTEC Inc.
- Do not soak the intraocular lens with any solution other than a sterile balanced salt or saline solution.
- Do not use the intraocular lens after the expiry date shown on the outside of the package. After this date, Lenstec cannot guarantee that the performance of the IOL will remain unchanged.
- Handle the intraocular lens carefully. Rough or excessive handling may damage the lens. A high level of surgical skill is required for intraocular lens implantation.
- A surgeon should have observed and/or assisted in numerous surgical implantations and successfully completed one or more courses on intraocular lenses, prior to attempting to implant lenses.
- Prior to implant, the surgeon must examine the intraocular lens to assure the correct haptic angle orientation, and the optic should vault anteriorly from the haptic footplates.
- All cases of lens removal must be reported to LENSTEC.
- Medical facilities utilizing this IOL, and its accessories (if any), must ensure proper disposal as medical waste.

## ADVERSE REACTIONS

The following adverse reactions have been reported following cataract extraction and implantation of an intraocular lens:

Endophthalmitis, hyphema, hypopyon, intraocular infection, lens dislocation, cystoid macular oedema, pupillary block, retinal detachment, secondary surgical intervention (excluding retinal detachment and posterior capsulotomy), corneal stromal oedema, iritis, cystoid macular oedema, raised intraocular pressure requiring treatment.

## HOW SUPPLIED/DIRECTIONS FOR USE

LENSTEC INC. PMMA intraocular lenses are ethylene oxide sterilized in a lens tray contained within two (2) sealed tyvek sterilizable peel pouches. The contents of the pouch are sterile unless the pouch has been opened or damaged. Remove the inner pouch from the outer pouch over or within a sterile field. Remove the lens tray from within the second pouch. Firmly hold the lens tray in the palm of hand and slide the lid open with your thumb; inspect the lens for any damage and debris. Grasp the lens by the haptics and rinse in a sterile balanced salt solution prior to implantation.

The lens box contains peelable labels that display the lens serial number, model name and number. These labels are designed to be affixed to the patient's hospital chart, physician's chart and the patient identification card. This card should be given to the patient as a permanent record of their implant.

## IMPLANTATION

The recommended surgical technique for implantation of Lenstec's PMMA intraocular lenses is with the use of a pair of forceps. Gently grasp the lens (by the optic-haptic junction) with the pair of forceps, and carefully introduce it into the anterior chamber (for LA-series IOLs) or the posterior chamber (for LS-series IOLs) of the eye. Gently adjust the IOL until the required position within the eye is obtained. Sutures are advisable if the incision section appears leaky.

In cases where a surgeon considers this technique to be impractical for a particular procedure, he/she is likely to be best advised to use an alternative implantation method in which his/her training and judgment dictate to give the best results to the patient.

## EXPLANATION

Explanation procedures may vary depending on patient condition and circumstances. The surgeon is therefore advised to use an explanation method which he/she determines will provide the most favourable patient outcomes.

## CALCULATION OF LENS POWER

It is recommended that the surgeon use a power calculation method in which he is most familiar and comfortable with. In general, the power of the lens for each patient can be estimated from prior refractive error or calculated from the corneal radius, depth of the anterior chamber and axial length of the eye according to formulas in corresponding literature.

*NOTE: The 'A' Constant and ACD values printed on the outside of the package are estimates only. It is recommended that the surgeon determine his/her own values based on their individual clinical experience.*

## ADVERSE EVENT REPORTING

Adverse reactions and/or potentially sight-threatening complications that may reasonably be regarded as lens related that were not previously expected in nature, severity, or incidence must be reported to LENSTEC.

## IMPLANT CARD

- All fields present on the ID card must be completed by the healthcare institution/provider. Please use a marker to complete empty fields.
- It is the responsibility of the healthcare institution/provider to attach the label sticker with the etched 'LENSTEC' logo to the reverse (unprinted) face of the patient ID card, and provide this to the patient as a record of their implant.

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