Table 3 - BEST CORRECTED DISTANCE VISUAL ACUITY at 1 Year (Form 5) Stratified by Age (Years)					
	< 60	60 to < 70	70 to < 80	≥ 80	
20/10 or better	0 / 36 (0%)	0 / 128 (0%)	0 / 155 (0%)	0 / 47 (0%)	
20/16 or better	6 / 36 (16.7%)	12 / 128 (9.4%)	4 / 155 (2.6%)	2 / 47 (4.3%)	
20/20 or better	26 / 36 (72.2%)	85 / 128 (66.4%)	78 / 155 (50.3%)	22 / 47 (46.8%)	
20/25 or better	33 / 36 (91.7%)	110 / 128 (85.9%)	121 / 155 (78.1%)	33 / 47 (70.2%)	
20/30 or better	36 / 36 (100%)	125 / 128 (97.7%)	143 / 155 (92.3%)	41 / 47 (87.2%)	
20/40 or better	36 / 36 (100%)	127 / 128 (99.2%)	152 / 155 (98.1%)	45 / 47 (95.7%)	
20/50 or better	36 / 36 (100%)	127 / 128 (99.2%)	153 / 155 (98.7%)	47 / 47 (100%)	
20/60 or better	36 / 36 (100%)	127 / 128 (99.2%)	154 / 155 (99.4%)	47 / 47 (100%)	
20/80 or better	36 / 36 (100%)	127 / 128 (99.2%)	154 / 155 (99.4%)	47 / 47 (100%)	
20/100 or better	36 / 36 (100%)	127 / 128 (99.2%)	155 / 155 (100%)	47 / 47 (100%)	
20/200 or better	36 / 36 (100%)	127 / 128 (99.2%)	155 / 155 (100%)	47 / 47 (100%)	
Worse than 20/200	0 / 36 (0%)	1 / 128 (0.8%)	0 / 155 (0%)	0 / 47 (0%)	
Not Reported	0	0	0	0	
Total	36	128	155	47	

Table 4 - Patient Population Softec HD™ PCIOL				Table 5			
n = 366 eyes in 366 study subjects with 1 year follow-up Cumulative Adverse Event Softec HD TM End PCIOL Grid			Persistent Adverse Event at 6 mths and/or 1 year	Softec HD™ PCIOL	FDA PCIOL Grid		
through 1 year	PCIOL Incidence	n = 300		Corneal Stromal Edema	0%	1.8%	
Cystoid Macular Edema	0.8%*	6.0%		Cystoid Macular Edema	0.8%*	2.2%	
Hypopyon	0%	1.8%		ojular madalar Edoma	0.0%	1.000	
Endophthalmitis	0%	1.0%	1	Intis Paised IOP Paguiring	0.3%	1.8%	
Dislocated Lens (from Posterior Chamber)	0%	1.0%		Treatment	0.3%	1.8%	
Pupillary Block	0%	1.0%	1	*Identical cases reported in persistent & cumulative CM **All unrelated to Softec HD™ PCIOL			
Retinal Detachment	0%	1.8%	1				
Secondary Surgical Intervention**	0.8%	2.6%					

Cumulative adverse events consist of all adverse events (AEs) that occurred at any point in postoperative follow-up during the first year after Softer HDTM PCIOL surgery.

Table 4 present all cumulative adverse events through the 1 year visit (330-420 days); Table 5, all persistent adverse events at 6 months (120-180 days) and 1 year visits. The overall incidence of cumulative and persistent IOL Grid adverse events in the Softee HD™ PCIOL Study Group (n = 560) was 2.2% (CME 0.6% secondary surgical interventions 0.8%, initis 0.3% and raised OP requiring treatment 0.3%). Non-IOL Grid AEs included 9 haptic break AEs at the time of the initial surgery and 1 sub-retinal hemorrhage.

EXPIRATION DATE The expiration date on the lens package is the sterility expiration date. Do not use the IOL after the expiration date

RETURNS POLICY Contact your Lenstec representative regarding the return goods policy. Return the lens with full identification and the reason for the return. Label the return package as a biohazard.

Lade the Pittum package is a unitazitu. PATIENT REGISTRATION AND REPORTING A Patient Identification Card is included in the package. This is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye practitioner the patient consults in future. Self-adhesive lens identification labels are provided for use on the Patient Identification Card and other clinical records. Adverse events/compliants that may reasonably be regarded as lens-related and that were not previously expected in nature, severily, or degree of incidence should be reported to Lenstec (Rahados) (n., Arport Commercial, Pliqtim Raad, Christ Church, Barbados: Tel: 246-420-6795 • Fax: 246-420-6797; Email: <u>Feedback@Lenstec.com</u> or contact your Lenstec representative.

BIBLIOGRAPHY

- LIOCRAPHY Holladay, J. Quality of Vision: essential optics for the cataract and refractive surgeon. Pgs 27-35 Apple DJ. Arthur S, et al. Effect of heperin surface modification in reducing silicone oil adherence to various intraocular lenses. Journal of Cataract & Refractive Surger, October 2001, Vol 27 No 10, 1662-1669.
 Hayashi H, Hayashi K, et al. Qualitative comparison of posterior capsule opacification after polymethylmethacrylate, silicone and soft acrylic intraocular lens implantation. <u>Arch. Ophthalmol.</u> 1998; 116: 1579-82.
 Apple DJ, Fedman JL, et al. Interversible silicone oil adhersion to silicone intraocular lenses. A clinicopathologic analysis. <u>Ophthalmology</u> June 1996; Vol 103, 1555-1561.
 Pandey SK, Apple DJ, et al. Posterior Capsule Opacification: A Review of the Aetiopathogenesis, Experimental and Clinical Studies and Factors for Prevention. Indian J <u>Ophthalmol</u> 2004; 52:99-112.
 Apple DJ, Keinmann G, et al. A new classification of calcification of intraocular lenses. <u>Ophthalmology</u> Jan 2008, Volume 115, Issue 1, Pages 73-79.
 ISO 11979-7 Ophthalmic Implants Intraocular lenses Part 7: Clinical Investigations; 2006.

IMPLANT CARD

- All fields present on the ID card must be completed by the healthcare institution/provider. •
- It is the responsibility of the healthcare institution/provider to attach the label sticker with the etched 'LENSTEC' logo to the reverse (unprinted) face of the patient ID card, and provide this to the patient as a record of their implant.



河月 http://www.lenstec.com ISO13485 Registered Compa

LENSTEC (BARBADOS) INC., Airport Commercial Centre, Pilgrim Road, Christ Church BB17092, Barbados

INSTRUCTIONS FOR USE SOFTEC HD[™], SOFTEC HDO[™], SOFTEC I[™] AND SOFTEC HDM[™] POSTERIOR CHAMBER INTRAOCULAR LENSES (PCIOLs)







PI28 Rev 19





IMPORTANT NOTICE

It is highly recommended that the surgeon adheres to the recommendations, precautions, contraindications and warnings outlined in these instructions

CAUTION: Federal (U.S.) law restricts this device to the sale by or on the order of a physician.

CAUTION: Federal (U.S.) law restricts this device to the sale by or on the order of a physician. DEVICE DESCRIPTION The LENSTEC Softec series of Posterior Chamber Intraocular Lenses (PCIOLs) are ultraviolet absorbing, single-piece "C" loop intraocular lenses intended for the replacement of the human crystalline lens following phacoemulsification cataract removal. The Softec HD^{ID} possesses a balanced aspheric optic design (containing symmetrical aspheric anterior and posterior surfaces); Il is offered in the provide ranker of a posterior surfaces); Il is offered in the dispite optic aspheric optic design, with the +15.0 to +25.0 range being available in quarter diopter increments. The Softec HD^{ID}. The Softec II^{ID} has an optic designed with a balanced being ovoid rather than circular. It is offered in the dioptic power range of +5.0 to +36.0. The Softec HD^{ID}. The Softec II^{ID} has an optic designed with a balanced spherical optic. It is offered in the dioptic power range of +5.0 to +36.0. The Softec HD^{ID}. The Softec II^{ID} has an optic designed with a balanced spherical optic. It is offered in the dioptic power range of +5.0 to +36.0. The Softec HD^{ID}. The Softec ID^{ID}. Each of the Softer enses of intraorular lenses is manifectured from a mericial rande convolumer of Hutrohomitic Arvolic, with a nohumerisable II/I.

Tablete cannot upic intrasta and paper of the set of th

INDICATIONS FOR USE

The LENSTEC Softec HD™ Aspheric. Softec HDO™ Aspheric. Softec I™ and Softec HDM™ Posterior Chamber Intraocular Lenses (PCIOLs) are Intended for the replacement of the human crystalline lens following phaceemulsification cataract removal in adults over the age of 21. The lenses are indicated for capsular bag placement.

Cutside of general contraindications for ocular surgery, the following specific contraindications apply: Uncontrolled glaucoma, microphthalmia, chronic severe uveitis, retinal detachment, corneal decompensation, diabetic retinopathy, iris atrophy, perioperative complications, potentially foreseeable post-operative complications and other conditions which an ophthalmic surgeon might identify based on their experience.

WARNINGS

The implanting ophthalmic surgeon shall consider the following warnings, and identify a risk/benefit ratio prior to surgery: 1. Failure to follow the implantation instructions supplied with this lens could lead to mishandling and subsequent IOL damage prior to or during

implantation

- 2. 3. There is no clinical data to support placing this lens in the ciliary sulcus.
- 4. 5.
- 6. 7
- Impantation. There is no clinical data to support placing this lens in the cliary succes. Any posterior capsulotomy opening should be limited to approximately 4 mc. Consistent with other IOLs, there is an increased risk of lens dislocation and/or second support placing the capsular base is not intact of If there is significant transcesses. The OLS should not be implated by the capsular base is not intact of If there is significant transcesses. The effectiveness of ultraviolet light absorbing lenses in reducing the incidence of retinal disorders has not been established. As a precaution, patients should be informed that they should were surgisesses with UV protection when in swinght. The rate of cystol macular edema may increase with extracapsular bag placement of the haptics. Patients who present complications in the Contraindications' section of this document, non-ge related cataract, vitreous loss, ins atrophy, severe ensitexing, out less the norm-dage macular degeneration or suspected microbial limetion. Patients with any include, but is not limited to: persistent bleeding, significant its damage, uncontrolled positive pressure or significant vitreous prologes or loss. The implanting surgeon shall consider whether patients in who intraocular lens implantation would affect the ability to observe, diagnose or trad topositive segment diseases, should have any of the IOLs implanted. The implanting surgeon shall consider whether patients who have a distorted eve to previous trauma or developmental defects in which appropriate support of the IOL is no positive, should have any of the IOLs implanted. The implanting surgeon shall consider whether patients who have recurrent severe anterior or posterior segment inflammation or uveitis, should have any of the IOLs implanted.
- 8.
- 9.
- 10. 11
- 12. 13. 14.

- should have any of the IOLs implanted. Any circumstances which could lead to damage to the corneal endothelium during implantation should be avoided. Children under the age of 2 are not suitable candidates for intraocular lenses. Reuse of the IOL is shichly prohibited, as it raises serious safely and effectiveness concerns. LENSTEC does not provide cleaning/sterilization instructions. An improperty cleaned and/or sterilized IOL can cause significant damage to a patient's vision, due in part to rosso contamination induced infection. Once removed from its original packaging, the IOL can lose traceability. In the event an IOL is re-used, it is unlikely the user will know the correct expiry leads estain under or dioptic power. LENSTEC can not guarantee stability or proper function of either haptic or optic portions in the event that an IOL is re-used. Failure of either of these components can render the IOL ineffective.

PRECAUTIONS

- RECAUTIONS

 The IOL must be stored in dry conditions between 0°C (32°F) and 45°C (113°F).

 Do not attempt to re-use the lens. Do <u>not</u> autoclave or attempt to re-sterilize the lens. Lenses requiring re-sterilization should be returned to
 LENSTEC inc.
 Do <u>not</u> use the device if sterile packaging has been damaged or if there are traces of leakage on the bottle or pouch.
 Do <u>not</u> such the intraocular lens with any solution other than a sterile balanced sait solution or saline solution.
 Once packaging has been opened, the intraocular lens must be used immediately. The hydrophilic nature of the lens can cause the lens to
 absorb substances with which i comes into contact, such as disinfectants, medicines, lood cells, etc. This may cause a "Toxic Lens
 Syndrome". Rinse the lens carefully once removed from the glass vial.
 The lens must be immediated within 2 minutes following removal from its saline bath, as dehydration causes the lens material to

- The lens must be implanted within 2 minutes following removal from its saline bath, as dehydration causes the lens material to become brittle. The lens must be implanted in the capsular bag.

- The lens must be implanted using only injection systems validated for use with the IOLs. These include the Softac IOL Injection System (Lenster Inc), Viscojied 1.8th Injector Set, Model# (JP04450 (Medical AG) and the Softio injection system, Model# AS-9300 (ASICO). Do <u>not</u> use the intraocular lens after the expiration date shown on the outside package label. Handle the intraocular lens acterity, Rough handling or excessive handling may damage the lens. The surgeon must be aware of the risk of opacification of the intraocular lens, which may necessitate lens removal. NOTE: Although the LENSTEC hypothile intraocular lens has a salistactory history regarding lens opacification, there is a history of lens opacification with lenses from other manufacturers¹. The material used by LENSTEC, unlike the materials used by other manufacturers has not had any reported 'Adverse Events' due to material discloration, opacification andro other material tested deficiencies, which have caused opsolerative patient problems. Ophthalmic surgeons should keep in mind that there have been cases of reported opacification of hydrophilic IOLS, Most, if not all of these types of cases required explantation. All cases of lens removal must be reported to LENSTEC. **VEICENT**

HOW SUPPLIED

The IOLs are supplied in a 0.9% saline solution in a lens bottle contained within a sealed Tyvek sterilizable peel pouch and should only be opened under aseptic conditions

DIRECTIONS FOR USE

United index for Que a locate sterilized in a lens bottle contained within a sealed Tyvek startilizable peel pouch. The contents of the pouch/bottle are sterile unless the package is damaged or opened. NOTE: A blue IOL was used only to aid in contrasting the lens from the cartridge in the attached photograph.

INSTRUCTIONS FOR IMPLANTATION: SOFTEC PCIOL

Calculation of Lens Power

It is recommended that the surgeon use a lens power calculation method with which they are comfortable. In general, the power of the lens for each patient can be calculated from the keratometry measurements and axial length of the eye according to formulas in relevant literature. An A-Constant of 118.0 and an anterior chamber depth (ACD) of 5.10 should be used for the LENSTEC IOLs if an applanation A-Scan unit is used. This needs to be modified for other methods such as the IOL Master or Lenstar. Depending on the IOL power calculation formula being used by the physician, this value for use with the IOL Master or the Lenstar will change slightly. If using the SRK/T power calculation formula, this value should be 118.54 for the Softec HDO, 118.3 for the Softec HD and Softec I and 117.8 for the Softec HDM. Additional reference to this topic can be found at http://w

Pre-Surgical Preparation:

. er from IOI. Refractive Calculation Equation Holladay or SRK/T b. Determine the Expected Post-operative Target Refraction (SE).

SURGICAL TECHNIQUE Ensure capsulorhexis is up to 5.5 mm in diameter

b. Perform standard phacoemulsification technique

- b. Perform standard phacoemulatification technique.
 c. Unscrew the can from the glass viail. See Figure 1.
 d. Remove the Delivery System from the glass viail. See Figure 1.
 d. Remove the Delivery System from the glass viail. See Figure 1.
 e. Turn the Delivery System from the glass viail. See Figure 1.
 e. Turn the Delivery System from the glass viail. See Figure 1.
 f. Using tothless forceps, remove the lens by elimethe the applics or the optic, taking care not to cause damage to the lens. See Figure 4.
 g. If using the Medicei Viscolet 1.8¹⁷ Injector Set, Model# LP604350 or the ASICO Softip injection system. Model# AS-9300, proceed as directed in the respective injection systems. The ASICO Softip injection system. Model# AS-9300, proceed as follows: Prepare the injector cartidge (LENSTEC Softe OL Injection System) with viscolastic. Open the cartidge flaps and inject viscolastic down each side of the chamber, and into the tip (nosecone). See Figure 5.
 h. Laad the implant.

h. Load the implant.

- Obtain the Injector (LENSTEC Softec IOL Injection System) and make sure that the tip is exposed. Use the applicator to affix the silicone tip onto the injector tip and then pull its plunger back as far as it will go. See Figure 6. • Holding the flaps of the cartridge open as far as possible, place the lens in the cartridge chamber as indicated in Figure 7. Ensure that the

- The provide theprovide theprovide the provide theprovide the provide the prov
- will go
- Depress the injector plunger so that the silicone tip fits into the back of the cartridge chamber and advance it forward until you can just see the silicone tip in the cartridge tip (nosecone)

- The injector is now ready to use.
 Carefully introduce the loaded injector tip into the anterior chamber (bevel facing down to avoid touching the endothelium) until the opening of the cartidge is beyond the distal pupil margin. Gently inject the lens. Rotate the injector counterclockwise if necessary to ensure the IOL remains orientated correctly as it emerges from the cartidge. Ensure the leading haptic is in the bag and the lens haptic is orientated correctly. See Figure 11.
- correctly. See right? 11. k. Genity withdraw the cartridge from the eye as the trailing haptic emerges from the cartridge. I. Reconfirm that the anterior chamber is deep, and if not, introduce additional viscoelastic material. m. Using a tapered "pusher" insert the trailing haptic if protruding from the section and let it drop into the bag.
- n. Irrigate out the viscoelastic from the anterior chamber and from behind the IOL
- Hydrate the edges of the section to seal it. No sutures are normally required but if the section appears leaky or the chamber remains shallow. a suture may be advisable.

DETAILED DEVICE DESCRIPTION

Single Piece 26% Water Content HEMA (Hydroxyethyl methacrylate)



Index of refraction; 1.460

	Softec HD™, Softec I™, Softec HDM™	Softec HDO™
Optic Size	5.75mm (circle)	5.75 x 6.50mm (oval)
Optic Type	Equiconvex	Equiconvex
Length	12.00mm	12.50mm
Angulation	0 degrees	0 degrees
Construction	1 Piece	1 Piece
Position Holes	0 Holes	0 Holes
Optic/Haptic Material	HEMA (26% water content)	HEMA (26% water content)
A-Constant*	118.0**	118.0**

The lens specifications for the LENSTEC Softec series of Posterior Chamber Intraocular Lenses are as follows Intraocular Lens Specifications

*Guidelines for Calculation of Implant Powe **See above section titled 'Calculation Lens Power'. Diopter crements Diopter Softec HD™ Power Ranges Tolerance Applied* Softec I™ Power Ranges Tolerance Applied*** crement Offered In Offered In The Softec HDI™, Softec HDO™, Softec I™ and Softec HDM™ Posterior Chamber Intraocular Lenses are manufactured in the following dioptric ranges: +5.0 D to +10.0 D 1.0 D (± 0.25 D) +5.0 D to +10.0 D 1.0 D (± 0.25 D) +10 5 D to +14 5 D 05D (+ 0.25 D) +10.5 D to +30.0 D 05D (± 0.25 D) +15.0 D to +25.0 D 0.25 D (± 0.11 D) +31.0 D to +36.0 D 1.0 D (± 0.5 D) +25.5 D to +30.0 D 0.5 D (± 0.25 D) +31.0 D to +36.0 D 1.0 D (± 0.5 D) Diopter Diopter Softec HDO™ Tolerances Applied*** Softec HDM™ Tolerance crements ncrements Power Ranges Power Ranges Applied* Offered In Offered In +5.0 D to +10.0 D 1.0 D (± 0.25 D) +5.0 D to +10.0 D 1.0 D (± 0.25 D) +10.5 D to +14.5 D 0.5 D (± 0.25 D) +10.5 D to +14.5 D 0.5 D (± 0.25 D) ***Internal manufacturing/sorting tolerance +15.0D to +25.0 D 0.25 D (+ 0.11 D) +15.0D to +25.0D 0.25 D (± 0.11 D) +25.5D to +30.0 D 0.5 D (± 0.25 D) +25.5D to +30.0D 0.5 D (± 0.25 D) +31.0 D to +36.0 D 1.0 D (± 0.5 D) +31.0 D to +36.0 D 1.0 D (± 0.5 D)

COMPATIBILITY GUIDE

	IOL Injection Systems						
IOL Model		Softec		Viscoject	Softip		
	Validated for Use	Power range (D)	Validated for Use	Power range (D)	Validated for Use	Power range (D)	
Softec HD™	~	I-9011S/ LC1620: 5.0 to 26.0		5.0 to 26.0	~	E 0 to 26 0	
		I-9011S/ LC2420: 5.0 to 36.0	•			3.0 10 20.0	
Softec HDOTM		I-9011S/ LC1620: 5.0 to 20.0	~	N/A	×	NIA	
	•	I-9011S/ LC2420: 5.0 to 36.0	^			NVA	
Softec I™	>	I-9011S/ LC1620: 5.0 to 26.0		5.0 to 26.0		5.0 to 26.0	
		I-9011S/ LC2420: 5.0 to 36.0	•		•		
Softec HDM™	<	I-9011S/ LC16: 5.0 to 36.0	×	N/A	×	N/A	

CLINICAL OUTCOMES

CLINCAL OUTCOMES The multi-center U.S. Softec HD™ PCIOL Clinical Investigation was conducted at 8 clinical centers with Softec HD™ PCIOL implantations occurring between December 13. 2006 and June 9, 2006. One year postoperative follow-up provides documented evidence of the safety and effectiveness of the Softex HD™PCIOL for the indications for use stated in this physical nabeling.

Patient Population Three hundred and ninety eyes of 390 Three hundred and ninety eyes of 390 study subjects were implanted with the Softec HD[™] PCIOL. The Softec HD[™] Study Cohort consisted of 227 females and 163 males; 334 were Caucasian, 11 Black, 6 Asian, 4 Mixed and 35 "Other". The mean age for the study cohort was 20 a ware Caucas 70.8 years. One year follow-up was collected for 366 eyes of 366 study subjects

Table 1 - Patient Population - Softec HD™ PCIOL n = 390 eyes in 390 study subjects						
Patient Population	Population Description					
Mean Age (years)	70.8 yrs					
Patients with Pre-existing Macular	3.1%					
Other Patients with Pre-existing	30.5%					
Gender	Female Male	58.2% 41.8%				
Race	Caucasian Black Asian Mixed	85.6% 2.8% 1.5% 1.0%				

Visual Acuity Table 2 summarizes the postoperative visual acuity outcomes at the 1 year visit (330.420 days) for the Softee HDTM PCIOL Study Group who did not have a preoperative ocular pathology or postoperative macular degeneration ("Best Case" Cohort).

Table 3 for "All Eyes" Cohort in the Softec HD™ PCIOL Study Group. Note: 30 study subjects had YAG capsulatomies 12 months or earlier, 17 six months or less, YAG capsulatomy is anticipated to produce an improved BCVA outcome versus a pre-YAG outcome.

20/50 or better

20/60 or better 20/80 or better

20/100 or better

20/200 or better

Worse than 20/200

Not Reported

Total

32 / 32 (100%)

32 / 32 (100%)

32 / 32 (100%)

0/32(0%)

0

32

Other Patients with Pre-existing Conditions			30.5%		
Gender		Female Male	58.2% 41.8%		
Race		Caucasian Black Asian Mixed Other	85.6% 2.8% 1.5% 1.0% 9.0%		
Table 2 - BE	AL ACUITY at 1 Ye by Age (Years)	ar (Form 5)			
	< 60	60 to < 70	70 to < 80	≥ 80	
20/10 or better	0 / 32 (0%)	0 / 118 (0%)	0 / 135 (0%)	0 / 42 (0%)	
20/16 or better	5 / 32 (15.6%)	12 / 118 (10.2%)	4 / 135 (3%)	2 / 42 (4.8%)	
20/20 or better	24 / 32 (75%)	79 / 118 (66.9%)	68 / 135 (50.4%)	21 / 42 (50%)	
20/25 or better	30 / 32 (93.8%)	100 / 118 (84.7%)	108 / 135 (80%)	31 / 42 (73.8%)	
20/30 or better	32 / 32 (100%)	115 / 118 (97.5%)	127 / 135 (94.1%)	39 / 42 (92.9%)	
20/40 or better	32 / 32 (100%)	117 / 118 (99.2%)	132 / 135 (97.8%)	42 / 42 (100%)	

32 / 32 (100%) 117 / 118 (99,2%) 133 / 135 (98,5%) 42 / 42 (100%)

32 / 32 (100%) 117 / 118 (99.2%) 134 / 135 (99.3%) 42 / 42 (100%)

117 / 118 (99.2%) 134 / 135 (99.3%)

135 / 135 (100%)

135 / 135 (100%)

0 / 135 (0%)

0

135

117 / 118 (99.2%)

117 / 118 (99.2%)

1/118(0.8%)

118

42 / 42 (100%) 42 / 42 (100%)

42 / 42 (100%)

0 / 42 (0%)

0

42