<http://www.lenstec.com>
ISO 13485 Registered Company

LENSTEC (BARBADOS) INC., Airport Commercial Centre, Pilgrim Road, Christ Church BB17092, Barbados

+1 (246) 420-6795

+1 (246) 420-6797

lenstecbarbados@lenstec.com

UK SALES: LENSTEC (BARBADOS) INC., Lenstec House, Unit 8, Mariner Court, Calder Park, Wakefield, WF4 3FL, UK

+44 (0) 1924 382 678

+44 (0) 1924 850 454

lenstecuk@lenstec.com

WORLD HEADQUARTERS: LENSTEC, INC., 1765 Commerce Ave N, Saint Petersburg, Florida 33716, U.S.A.

+1 (727) 571-2272


+1 (727) 571-1792

lenstec@lenstec.com

DISPOSABLE INJECTION SYSTEM

Model	Tip Angle (°)	Tip Diameter (mm)
LC16	45	1.6
LC24	45	2.4
Cart45S	45	1.6

STERILE EO

R_{only}

DEVICE DESCRIPTION

Lenstec cartridges are sterile, single use, disposable accessories, used to fold the intraocular lens prior to implantation. The devices possess a "twin flap" design with a small tip to allow the surgeon to fold and insert the intraocular lens through a small incision.

INTENDED USE

The Lenstec family of cartridges is intended to aid qualified surgeons in the implantation of a foldable intraocular lenses into the eye following cataract extraction. These products are required to successfully facilitate implantation through a small incision without damaging the intraocular lens.

IMPORTANT NOTICE

It is highly recommended that the surgeon adhere to the recommendations, contraindications and warnings outlined in these instructions.

RISKS

The potential risks of using the Lenstec Disposable Injection Systems are as follows: Capsular damage, Endophthalmitis, Toxic anterior segment syndrome (TASS), Uveitis, Iritis, Anterior capsule fibrosis

WARNINGS

- Cartridge injection systems manufactured by Lenstec are intended to be used only with IOLs manufactured by Lenstec. A **Compatibility Guide** is provided at the end of this document.
- Surgical instructions for use of the cartridge injection system with Lenstec's IOLs are provided in the IOL instructions for use.
- Do not use the cartridge after the expiration date shown on the outside package label.
- Proper procedure is the responsibility of the individual surgeon.
- Medical facilities utilizing these devices, and their accessories (if any), must ensure proper disposal as medical waste.

A, C, D

The titanium injector and lens loader are reusable (following decontamination and sterilization). The injector and lens loader are supplied unsterile and must be cleaned and sterilized prior to initial use.

Manual cleaning prior to and after use:

- Rinse injector or lens loader with deionized water and then clean the device and its crevices using a soft surgical/laboratory brush moistened with alcohol¹. Care should be taken to avoid damage to the extreme tip.
- Inspect the instrument under magnification (10x) to check for gross contamination, wear or damage. Do not use/reuse if wear or damage are apparent. Re-clean if gross contamination is visible¹.

Sterilization method(s) prior to and after use:

- Place unwrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 135°C for an exposure time of 3 minutes and minimum dry time of 0-1 minute.
- OR
- Place wrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 132°C for an exposure time of 15 minutes and dry time of 15-30 minutes.
- OR
- Place wrapped instruments onto trays and steam sterilize using gravity displacement sterilizer for a minimum 121°C for an exposure time of 30 minutes and dry time of 15-30 minutes.

B

- Re-use of the disposable injector is strictly prohibited, as it raises serious safety and efficacy concerns.
- LENSTEC does not provide cleaning/sterilization instructions. An improperly cleaned and/or sterilized disposable injector can cause significant damage to a patient's vision and health, due in part to cross contamination induced infection.
- Once removed from its original packaging, the disposable injector can lose traceability. In the event the disposable injector is re-used, it is unlikely the user will know the correct expiry date or serial number.
- LENSTEC cannot guarantee stability or proper function of the disposable injector in the event that an injector is re-used.

QUALITATIVE & QUANTITATIVE DATA FOR THE CARTRIDGE INJECTION SYSTEMS

Cartridge Injection System: Polypropylene cartridge. The devices have been tested and proven safe in accordance with ISO 10993-5, ISO 10993-7, ISO 10993-10, ISO 10993-11.

EXPIRATION DATE

The expiration date on the lens package is the sterility expiration date.

RETURNS POLICY

Contact your Lenstec representative regarding the return goods policy.

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

Bibliography

1. Chang, D. F., Mamalis, N., "Guidelines for the Cleaning and Sterilization of Intraocular Surgical Instruments", Ophthalmic Instrument Cleaning & Sterilization Guidelines (OICS) Task Force

COMPATIBILITY GUIDE

MODEL	POWERS (D)	ACCESSORIES
 Softec HP1™	+10.0D to +30.0D (LC16; Cart45S)	
 Softec I™ Softec HD™ Softec HDY™	Softec I -5.0D to +26.0D (LC16; Cart45S) -5.0D to +36.0D (LC24) Softec HD, HDY +5.0D to +26.0D (LC16; Cart45S) +5.0D to +36.0D (LC24)	
 Softec HDO™	+5.0D to +26.0D (LC16; Cart45S) +5.0D to +36.0D (LC24)	
 Softec HDM™	+5.0D to +36.0D (LC16)	
 SBL-3™ SBL-2™	+10.0D to +26.0D (LC16; Cart45S) +10.0D to +36.0D (LC24)	
 Tetraflex HD™	+5.0D to +26.0D (LC16; Cart45S) +5.0D to +36.0D (LC24)	

<u>I9011S Reusable Injector</u>	<u>Disposable Injectors</u>
LC16 (1.6mm)	LC16I (1.6mm)
LC24 (2.4mm)	LC24I (2.4mm)
Cart45S (1.6mm)	Cart45SI (1.6mm)