


Pre-Loaded Injector



1/4
DIOPTER


Pre-Loaded Injection System:
All in One!

- » Reduces handling and loading errors
- » Minimizes the chance of infection
- » Saves time - quicker operating room turnover

The TetraflexHD and PLI V2 are not currently for sale in the US.




Pre-Loaded Injector



1 Place the PLI Tip into BSS or saline, then insert the tip into the slot



2 Remove the red button from the PLI



3 Firmly depress the clear tab



4 Advance the lens to the tip using the plunger and ensure on all sides that no part of the lens is trapped under the cushion

World Headquarters
Lenstec, Inc.
1765 Commerce Ave. N.
St. Petersburg, FL 33716
USA

Tel: 727-571-2272
Fax: 727-571-1792
Email: lenstec@lenstec.com

Lenstec Barbados
Lenstec Barbados Inc.
Airport Commercial Centre
Pilgrim Road, Christ Church
BARBADOS BB17092

Tel: 246-420-6795
Fax: 246-420-6797
Email: lenstecbarbados@lenstec.com

Lenstec UK
Lenstec Barbados Inc.
Lenstec House
Unit 8, Mariner Court
Calder Park, Wakefield
West Yorkshire WF4 3FL
England

Tel: +44 (0)1924 382 678
Fax: +44 (0)1924 850 454
Email: lenstecuk@lenstec.com



www.lenstec.com