


# Softtec HD3



**Bringing a  
Clear View of  
the World**



## The Lenstec Softtec HD3 IOL:

- » Three-piece lens
- » Patented bi-aspheric design
- » Equiconvex optic (+1.0D to +36.0D), Meniscus (-5.0D to 0.0D)
- » Fully-angulated lens gives it more reliable posterior position

# Softec HD3

## Softec HD3 Technical Specifications†

<b>Optic Size</b>	6.00 mm
<b>Optic Type</b>	Bi-Aspheric, Equiconvex(+1.0 to +36.0D) Meniscus (-5.0 to 0.0D)
<b>Length</b>	13.25 mm
<b>Haptic Style</b>	Modified C
<b>Angulation</b>	5 Degrees
<b>Positioning Holes</b>	0
<b>Construction</b>	3 Pieces
<b>Optic Material</b>	Acrylic Hydrophilic

## Constants (Optical Interferometry):\*


<b>Immersion</b>	A= 118.13
<b>SRK/T</b>	A = 118.13
<b>Holladay1</b>	sf = 1.30
<b>Hoffer Q</b>	pACD = 5.04

\*i.e. Using IOL Master, LENSTAR Optical Biometers

## Diopter Range

<b>Whole:</b>	-5.0 to +36.0
<b>Half:</b>	+10.5 to +29.5

†The A-Constant and ACD values printed are estimates only. It is recommended that the surgeon determine his or her own values based on their individual clinical experience.



### World Headquarters

Lenstec, Inc.  
1765 Commerce Ave. N.  
St. Petersburg, FL 33716  
USA

Tel: 727-571-2272  
Fax: 727-571-1792  
Email: lenstec@lenstec.com

### Lenstec Barbados

Lenstec Barbados Inc.  
Airport Commercial Centre  
Pilgrim Road, Christ Church  
BARBADOS BB17092

Tel: 246-420-6795  
Fax: 246-420-6797  
Email: lenstecbarbados@lenstec.com

### Lenstec UK

Lenstec Barbados Inc.  
Lenstec House  
Unit 8, Mariner Court  
Calder Park, Wakefield  
West Yorkshire WF4 3FL  
England

Tel: +44 (0)1924 382 678  
Fax: +44 (0)1924 850 454  
Email: lenstecuk@lenstec.com



**LENSTEC**  
www.lenstec.com