



Insertion Technology



I9011S Injector
Reusable



Cartridge: LC 16 Series
Bevel: 45 or 20 degrees



**LC16I or 1620I Disposable
Injector/Cartridge**



**Lens
Loader**



**Lens
Loader**



Cartridge: LC 16 Series
Bevel: 45 or 20 degrees



I-9012FS Injector
Reusable



Cartridge: LC 16 Series
Bevel: 45 or 20 degrees



Note: The Lens Loader is to be used with all injectors.



www.lenstec.com



**LENS
TEC**

Insertion Technology



LC I Series Disposable Injection System *In a Convenient 10-Pack*

A dispenser for single use injectors. Safe, convenient and cost effective solution to micro-incision simplicity. Increases safety, reduces preparation time and inventory cost. Maximizes control.



LensVisc: Viscoelastic Solutions from Lenstec

Manufactured to high medical standards, LensVisc solutions are packaged in convenient, sterile, single-use syringes. LensVisc solutions are biocompatible and completely transparent with no

LensVisc HA Viscoelastic Solution

- 1.0 ml of HA (Hyaluronic Acid 1.4% or 2%)
- 27 gauge cannula
- LensVisc HA is a totally synthetic material and contains no animal by-products



Not for sale in U.S.

air bubbles and offer ideal viscosity for ophthalmic surgery, protecting the capsular bag from collapse as well as reducing eye dehydration and inflammation.

LensVisc MC Viscoelastic Solution

- 1.5 ml of MC (Hydroxypropyl Methylcellulose 2%)
- 23 gauge cannula

Lenstec UK
 Lenstec Barbados Inc.
 Lenstec House
 Unit 8, Mariner Court
 Calder Park
 Wakefield WF4 3FL
 UK
 Tel: +44 (0)1924 382 678
 Fax: +44 (0)1924 850 454
 Email: lenstecuk@lenstec.com

Lenstec Barbados
 Lenstec Barbados Inc.
 Airport Commercial Centre
 Pilgrim Road
 Christ Church BB17092
 Barbados
 Tel: 246-420-6795
 Fax: 246-420-6797
 Email: lenstecbarbados@lenstec.com

World Headquarters
 Lenstec Inc.
 1765 Commerce Ave. N.
 St. Petersburg, FL 33716
 USA
 Tel: 727-571-2272
 Fax: 727-571-1792
 Email: lenstec@lenstec.com