

### Lens Specification

|                           |                                                                      |
|---------------------------|----------------------------------------------------------------------|
| <b>Optic Size:</b>        | 6.0mm                                                                |
| <b>Optic Type:</b>        | Bi-aspheric,<br>Equicovex (+1.0 to +36.0D) - Meniscus (-5.0 to 0.0D) |
| <b>Length:</b>            | 13.25mm                                                              |
| <b>Haptic Style:</b>      | Modified C                                                           |
| <b>Angulation:</b>        | 5 Degrees                                                            |
| <b>Positioning Holes:</b> | 0                                                                    |
| <b>Construction:</b>      | 3 Piece                                                              |
| <b>Optic Material:</b>    | Acrylic Hydrophilic (26% Water Content)                              |

### A-Constant / Biometry

|                   |             |
|-------------------|-------------|
| <b>Immersion:</b> | A= 118.13   |
| <b>SRK/T</b>      | A= 118.13   |
| <b>Holladay1:</b> | sf = 1.30   |
| <b>Hoffer Q:</b>  | pACD = 5.04 |

### Power Range

|                             |                |
|-----------------------------|----------------|
| <b>Whole Power Range:</b>   | +5.0 to +36.0  |
| <b>Half Power Range:</b>    | +10.0 to +30.0 |
| <b>Quarter Power Range:</b> | +15.0 to +25.0 |

**Please note: Available power ranges may vary on location and not all products are approved in all countries. Please contact us for availability.**

A constant and A/C Depth figures shown are strictly guidelines for the calculation of implant power. Lenstec recommends that the surgeons develop their own values based on technique, measuring equipment and desired postoperative results.